




DONJOY®

DonJoy ActyLight®

Fit-and-forget ankle support

enovis™





DONJOY®

DonJoy ActyLight®

Fit-and-forget ankle support

The convenience and comfort of the DonJoy ActyLight® ankle support means patients will forget they're even wearing it. Its step-in design and lockable lace system provide a fast, secure fit that's complemented by a pair of removable anatomical stays. Their simple adaptability enables both acute and chronic use, so ActyLight can support patients throughout their treatment.

Intended use/Indications

DonJoy ActyLight is designed to provide support and stabilization to the ankle.

- Initial management of chronic ankle instabilities
- Mild to moderate sprains of the lateral ligament
- Isolated syndesmotic injuries (Grade 1)
- Return to sports, prevention of re-injury
- Ligament instability
- Supination prevention


Features and benefits



Step-in design makes donning and doffing quick and easy



Secure lace-locking mechanism enables precise tightening, helps maintain a consistent fit, and allows quick release



Lightweight, breathable fabrics help the foot stay cool and dry to aid compliance



Figure-8 straps provide inversion/eversion control



Anatomically contoured bilateral stays provide additional stability and protection and can be repositioned to accommodate the malleoli




Fits inside a trainer comfortably allowing patients to stay active



Modular design means the support can be adapted easily for acute or chronic use by removing the stays







How it Works

ActyLight's step-in-sleeve design provides proprioceptive feedback and works with supportive elements to stabilize the ankle and restrict excessive range of movement.

- **Bilateral semi-rigid supports** provide mechanical stability to prevent ankle inversion/eversion
- **Figure-8 straps** provide further ankle inversion/eversion control
- **Forefoot mesh** made from breathable material to provide proprioceptive and compressive support for the foot



ORDERING INFORMATION

PART NO.	HEEL CIRCUMFERENCE	EU SHOE SIZE	SIZE
82-0054-1	≤ 27 cm	32 - 36	XS
82-0054-2	27 - 30 cm	36 - 39	S
82-0054-3	30 - 34 cm	39 - 43	M
82-0054-4	≥ 34 cm	43 - 47	L



**Asia Pacific**

1905, Tower 2, Grand Central Plaza
Shatin
HONG KONG
T: 852 31051415
E: marketing.ap@enovis.com

Australia – New Zealand

PO Box 6057
Frenchs Forest DC
NSW 2086
AUSTRALIA
T: 1300 66 77 30
E: customerservice.au@enovis.com

Benelux

Kleinhoefstraat 5, bus 39
Geel 2440
BELGIUM
T: +32 (0)14 248090
E: orders.bx@enovis.com

Canada

6485 Kennedy Rd
Mississauga
Ontario, L5T2W4
CANADA
T: 1-866-866-5031
E: canada.orders@enovis.com

China

A312 SOHO ZhongShan Plaza
1055 W. ZhongShan Rd
Shanghai, 200051
CHINA
T: (8621) 6031 9989
E: shujun.tang@enovis.com

France

3 rue de Bethar
Centre Européen de Frêt
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
E: sce.cial@enovis.com

Germany

Bötzingen Straße 90
79111 Freiburg
GERMANY
T: +49 (0) 180 1 676 333
E: kundenservice@enovis.com

Italy

Via Leonardo da Vinci 97
Trezzano sul Naviglio
Milan 20090
ITALY
T: +39 02 835 98001
E: vendite@enovis.com

Nordic

Murmansgatan 126
212 25 Malmö
SWEDEN
T: +46 40 39 40 00
E: info.nordic@enovis.com

South Africa

Unit 6, 2 on London, 2 London Circle
BrackenGate Business Park
Brackenfell, 7560
SOUTH AFRICA
T: +27 (0) 21 276 2968
E: info.southafrica@enovis.com

Spain

Av.Cornella 144. 1-4
08950 Esplugues de Llobregat
Barcelona
SPAIN
T: +34 93 480 32 02
E: ventas@enovis.com

UK & Ireland

1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UK
T: +44 (0)1483 459 659
E: ukorders@enovis.com

USA

2900 Lake Vista Drive
Lewisville
TX 75067
U.S.A.
T: 760 727 1280
E: customercare@enovis.com



Creating Better Together™

enovis-medtech.eu

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-000-00159475-EN - Rev A

Copyright © 2023 DJO, LLC, a subsidiary of Enovis Corporation