

DONJOY®

DonJoy Spinostrap

BACK SUPPORT, IN HAND



enovis™

DONJOY®

DonJoy Spinostrap

Back support, in hand.

DonJoy Spinostrap puts prescribers and patients in control of back support. Its adjustable shoulder straps, bendable exoskeletal dorsal frame, and customizable lumbar belt enable comfortable individualized support and relief of the thoracic spine. Thanks to this semi-rigid back brace's range of intuitive features, back support and kyphotic postural correction are well within reach.

INTENDED USE/INDICATIONS

DonJoy Spinostrap is designed to provide support and stabilization and partial immobilization (restriction of movement) of the spine.

- Support of stable spinal fractures
- Relief of pain and facilitation of correction of postural dysfunction from dorsal hyper-kyphosis
- Postural dysfunction associated with conditions such as Parkinsons disease
- Straightening and support of the spine in the case of muscular insufficiency/weakness
- Pain relief and support of spinal osteoarthritis



FEATURES AND BENEFITS



Semi-rigid exoskeletal dorsal frame: to provide support and aid relief of the thoracic spine, aluminum core can be molded to fit the individual curvature of the back



Repositionable strap guide: can be placed at different heights on the dorsal frame to help control the corrective action on kyphotic posture



Adjustable lumbar belt: uses accessible strap tabs and finger pockets to provide adjustable support for the spine while leaving the abdomen free. Can be cut to required length



Adjustable shoulder straps: use fixable hook-and-loop strap tabs to provide adjustable tension, aiding correction of kyphotic posture by pulling the shoulders back and straightening the back



Crocodile strap closures: allow customization of shoulder strap and lumbar belt length for a convenient individualized fit



Axillary shoulder pads: shaped to be flat and anatomical, they can be positioned to increase the wearing comfort for the underarm area



Adjustable sternal strap: can be attached to provide additional comfort in the axillary area



Slim and low-profile design: fits under clothing

Lightweight and breathable: helps ensure patient comfort
Hypoallergenic and soft materials: to help avoid skin irritation. PVC and Latex free

ORDERING INFORMATION

PART NUMBER	SPINE LENGTH C7 – S1		SIZE
	CM	INCHES	
82-0057-1	42 – 52 cm	16½ – 20½"	S
82-0057-2	52 – 62 cm	20½ - 24½"	M
82-0057-3	62 – 72 cm	24½ - 28½"	L

**Asia Pacific**

1905, Tower 2, Grand Central Plaza
Shatin
HONG KONG
T: 852 31051415
E: marketing.ap@enovis.com

Australia – New Zealand

PO Box 6057
Frenchs Forest DC
NSW 2086
AUSTRALIA
T: 1300 66 77 30
E: customerservice.au@enovis.com

Benelux

Kleinhoefstraat 5, bus 39
Geel 2440
BELGIUM
T: +32 (0)14 248090
E: orders.bx@enovis.com

Canada

6485 Kennedy Rd
Mississauga
Ontario, L5T2W4
CANADA
T: 1-866-866-5031
E: canada.orders@enovis.com

China

A312 SOHO ZhongShan Plaza
1055 W. ZhongShan Rd
Shanghai, 200051
CHINA
T: (8621) 6031 9989
E: shujun.tang@enovis.com

France

3 rue de Bethar
Centre Européen de Frêt
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
E: sce.cial@enovis.com

Germany

Bötzingen Straße 90
79111 Freiburg
GERMANY
T: +49 (0) 180 1 676 333
E: kundenservice@enovis.com

Italy

Via Leonardo da Vinci 97
Trezzano sul Naviglio
Milan 20090
ITALY
T: +39 02 835 98001
E: vendite@enovis.com

Nordic

Murmansgatan 126
212 25 Malmö
SWEDEN
T: +46 40 39 40 00
E: info.nordic@enovis.com

South Africa

Unit 6, 2 on London, 2 London Circle
BrackenGate Business Park
Brackenfell, 7560
SOUTH AFRICA
T: +27 (0) 21 276 2968
E: info.southafrica@enovis.com

Spain

Av.Cornella 144. 1-4
08950 Esplugues de Llobregat
Barcelona
SPAIN
T: +34 93 480 32 02
E: ventas@enovis.com

UK & Ireland

1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UK
T: +44 (0)1483 459 659
E: ukorders@enovis.com

USA

2900 Lake Vista Drive
Lewisville
TX 75067
U.S.A.
T: 760 727 1280
E: customercare@enovis.com

enovis™*Creating Better Together™***enovis-medtech.eu**

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-000-00171850-EN - Rev A

Copyright © 2023 DJO, LLC, a subsidiary of Enovis Corporation